



Shipment of samples:

Register the samples electronically in PatoLink. *Contact us if you have questions about access or use of PatoLink.*

If the entire kit is not used, it can be stored in room temperature so that the remaining tubes can be used at a later point (note: expiry date).

Send the samples in a padded envelope with an ice pack, by overnight express to PatoGen's laboratory.

Use the link below to order a return label for sample submission:

[Sending samples - PatoGen](#)

NO: PATOGEN AS - Rasmus Rønnebergs gt. 21, 6002 Ålesund

UK: PATOGEN - Malin House - European Marine Science Park - Dunbeg - PA37 1SZ Oban